

MVIEWM.005A



# 71308  
10/13/01  
PATENT

IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

|           |   |  |                     |
|-----------|---|--|---------------------|
| Applicant | : | Sherman, et al.  | Group Art Unit 1652 |
| Appl. No. | : | 09/501,730   |                     |
| Filed     | : | February 10, 2000  |                     |
| For       | : | <b>AGGREGATE-FREE URATE<br/>OXIDASE FOR<br/>PREPARATION OF NON-<br/>IMMUNOGENIC POLYMER<br/>CONJUGATES</b> |                     |
| Examiner  | : | Pak, Y.  |                     |

**DECLARATION OF L. DAVID WILLIAMS, Ph.D.**  
**SUBMITTED UNDER 37 C.F.R. § 1.131**

Assistant Commissioner for Patents  
Washington, D.C. 20231

Dear Sir:

1. I, L. David Williams, Ph.D., am a named inventor on the above-captioned patent application.
2. I have read and understood the specification and the claims of the above-captioned patent application.
3. I have read and understood Caliceti et al. *Bioconjugate Chem.* **1999**, *10*, 638-646, which has a web publication date of June 2, 1999.
4. Exhibit A, attached hereto, is a copy of several pages of one of my laboratory notebooks.
5. Pages 1 and 2 of Exhibit A show the results of a purification procedure that I conducted prior to June 2, 1999. The tracings shown on pages 1 and 2 are ion-exchange and size-exclusion chromatograms, respectively, of the product of my purification steps. The areas

Appl. No. : 09/501,730  
Filed : February 10, 2000

under the peaks of the chromatograms on pages 1 and 2 show the relative abundance of several aggregates of uricase.

6. The areas under the peaks assigned to tetramers are much larger than the area under the peaks assigned to octamers or other aggregates. These data demonstrate the purification of uricase that is substantially free of aggregates larger than tetramers.
7. Therefore, prior to June 2, 1999 I had invented and was in possession the subject matter claimed in the above-captioned patent application.
8. I declare that all statements made herein are true, and that all statements made upon information and belief are believed to be true, and further, that these statements were made with the knowledge that willful, false statements and the like so made are punishable by fine or imprisonment, or both, under 18 U.S.C. § 1001, and that willful, false statements may jeopardize the validity of the application, or any patent issuing thereon.

October 4, 2001

Date

L. David Williams

L. David Williams, Ph.D.

S:\DOCS\SKT\SKT-1775.DOC  
100301

linear5\_M1\_UV\_280nm\_01 — linear5\_ConcB\_01 — fractions\_1

19-52

AU

Edit - Run Description

## METHOD

Method base : ml

\*\*\*\*\* Main method \*\*\*\*\*

|                     |       |    |
|---------------------|-------|----|
| 0.00 CONC_B         | 0.0   | () |
| 0.00 FLOW           | 1.0   | () |
| 0.00 RECORDER_SPEED | 0.50  | () |
| 130.00 CONC_B       | 0.0   | () |
| 250.00 CONC_B       | 100.0 | () |
| 270.00 CONC_B       | 100.0 | () |

Sample DTG-Ply-KS-dt, 80%, 0.2u filtered

Loop load thru pump

Column MuxQ H210/10

Eluent A 0.1M pH 10.3 NaCO3

Eluent B 0.6M NaCl in Buffer A

Gradient Linear 0-100% C from 130 to 250 ml

Flow rate 1 mL/min

Fraction size 4 mL

## RUN LOG

Method file name LINEARS.FMT

Run by LDW 8:52:26 am

Config: Method base is ml

Config: Injection valve is valve No. 1

Config: Pump A/B is P-500, calibration 110.00

Config: Pump C is P-1, calibration 150.00

Config: UV-M -&gt; Monitor 1, Input range = 100 mV

Method start. Result name LINEARS.FRS

|        |         |            |             |
|--------|---------|------------|-------------|
| 0.00   | Auto:   | Tube No 1  |             |
| 0.00   | Auto:   | Pause      | 10:13:02 am |
| 80.50  | Auto:   | Continue   | 10:16:54 am |
| 130.02 | Auto:   | Tube No 2  |             |
| 134.02 | Auto:   | Tube No 3  |             |
| 138.01 | Auto:   | Tube No 4  |             |
| 142.02 | Auto:   | Tube No 5  |             |
| 146.01 | Auto:   | Tube No 6  |             |
| 150.01 | Auto:   | Tube No 7  |             |
| 154.01 | Auto:   | Tube No 8  |             |
| 157.98 | Auto:   | Tube No 9  |             |
| 161.99 | Auto:   | Tube No 10 |             |
| 165.99 | Auto:   | Tube No 11 |             |
| 169.97 | Auto:   | Tube No 12 |             |
| 170.05 | Manual: | Hold       | 11:46:34 am |
| 173.96 | Auto:   | Tube No 13 |             |
| 177.96 | Auto:   | Tube No 14 |             |
| 181.95 | Auto:   | Tube No 15 |             |
| 182.03 | Manual: | Continue   | 11:58:33 am |
| 185.96 | Auto:   | Tube No 16 |             |
| 189.96 | Auto:   | Tube No 17 |             |
| 191.97 | Manual: | Hold       | 12:08:30 pm |
| 193.95 | Auto:   | Tube No 18 |             |
| 197.94 | Auto:   | Tube No 19 |             |
| 201.94 | Auto:   | Tube No 20 |             |
| 205.94 | Auto:   | Tube No 21 |             |
| 208.00 | Manual: | Continue   | 12:24:32 pm |
| 209.94 | Auto:   | Tube No 22 |             |
| 213.93 | Auto:   | Tube No 23 |             |
| 217.92 | Auto:   | Tube No 24 |             |
| 221.93 | Auto:   | Tube No 25 |             |
| 225.93 | Auto:   | Tube No 26 |             |
| 229.92 | Auto:   | Tube No 27 |             |
| 233.91 | Auto:   | Tube No 28 |             |
| 237.91 | Auto:   | Tube No 29 |             |
| 241.90 | Auto:   | Tube No 30 |             |
| 245.90 | Auto:   | Tube No 31 |             |
| 249.89 | Auto:   | Tube No 32 |             |
| 253.89 | Auto:   | Tube No 33 |             |
| 257.87 | Auto:   | Tube No 34 |             |
| 261.86 | Auto:   | Tube No 35 |             |
| 265.87 | Auto:   | Tube No 36 |             |
| 269.87 | Auto:   | Tube No 37 |             |
| 273.85 | Auto:   | Tube No 38 |             |
| 277.85 | Auto:   | Tube No 39 |             |
| 281.84 | Auto:   | Tube No 40 |             |
| 285.84 | Auto:   | Tube No 41 |             |
| 298.02 | Auto:   | Method end | 1:54:38 pm  |

1.5

1.0

0.5

0.0

0

50

100

12

13

14

15

16

17

18

19

20

21

22

23

24

25

26

27

28

29

30

31

32

33

34

35

36

37

38

39

40

41

42

43

44

45

46

47

48

49

50

51

52

53

54

55

56

57

58

59

60

61

62

63

64

65

66

67

68

69

70

tetramers<sub>us</sub>octamer<sub>us</sub>aggregates<sub>us</sub>pool<sub>us</sub>

Continued on Page 53

EXHIBIT A PAGE 1 OF 4

Read and Understood By

W. H. Z. J. J.

Signed

Date

Signed

Date

Flush Mono Q HR10/10 column in FPLC system with 0.6 M NaOH in 90% NaCl, then water, then a blank run.

Buffer A = 0.1 M pH 10.3 sodium carbonate, Buffer B = 0.6 M NaCl / buffer A.

Start loading 0.2  $\mu$  filtered BTG-Pig-KS-dN, lot AO-58-034 rec'd [redacted], from water-washed Nalgene filter unit, load 90 mL @ 1 mL/min

gross filter wt 159.6 g

net wet wt -62.4 g

97.2 g

1 mL/min

1.1 mL water/fill

99.3 g

directly then Pump A. About 8.5 mL remaining

8:39:50 start load

8:44:14 pause at end of 1<sup>st</sup> piston cycle.

8:52:04 start Method LINEAR 5

10:12:38 pause method, stop load, start Buffer A

10:16:30 @ 80.5 mL continue

10:56:05 @ 120 mL


11:06:08 Start fraction collector, 4 min/fraction = 4 mL/fraction, at #2.

Hold gradient for three fractions at ~~41.6%~~ 33.3% B = 0.2 M NaCl

Hold gradient for 4 fractions at 41.6% B = 0.25 M NaCl

13:34:14 reach 100% B, 2:50 mL

13:54:14 bump off, end @ 270 mL



Superdex 200 HR10/30

Buffer: 10 mM Carbonate + 100 mM NaCl, pH 10.1; 0.5 mL/min

UV2000 B 214nm = UV276nm

Continued on Page 54

EXHIBIT A PAGE 2 OF 4

Read and Understood By

Mark Sifer

Signed

Date

Signed

Date

Remove 180  $\mu$ l returns of 19-53-Fract's 17-22 into HPLC vials w inserts. Also prepare dilutions as follows

|    |      |                                 |
|----|------|---------------------------------|
| 17 | 1x   | 360 $\mu$ l return              |
| 18 | 1/8  | 50 $\mu$ l + 350 $\mu$ l buffer |
| 19 | 1/12 | 30 $\mu$ l + 330 $\mu$ l water  |
| 20 | 1/8  | 50 $\mu$ l + 350 $\mu$ l water  |
| 21 | 1/5  | 70 $\mu$ l + 280 $\mu$ l water  |
| 22 | 1/3  | 100 $\mu$ l + 200 $\mu$ l water |

Pool remaining 17-22 after return using fract 17 as wash of other to Pool wt = 21.565g

Remove 2 mL Pool return  
 $(21.565 - 2) \times 10.6 \frac{\text{mg}}{\text{mL}} = 20.7 \text{ mg target}$

Naigh out 233 mg NPC-10K-PEG (Sbeamwater LG-128-01) and add 0.932 mL 1 mM HCE at 18:16. Momentary yellow color on solids before they are well suspended. Residual base in NPC-PEG? Vt of added HCE = 0.941g. Add 1.099g PEG soln dropwise with stirring to the 19.5 mL Pool, ice bath, over 1-2 min. 1 mL PEG soln weighed 1.028g. Calcd NPC PEG added =  $233 \text{ mg} \times \frac{1.099}{0.941 + 233} = 21.8 \text{ mg}$   
 Covered + refrigerated

Remove  $\mu$ l sample of rxn for HPLC

Prepared 10K Midbee ultrafilter

by rinsing with water, then 0.1 M NaOH

Rinsed with water, then 0.1 M pH 10.3 buffer (left in lines, i.e. ~3 mL) and began ultrafiltering round 1 reaction at 16:10. By 16:23 of about 5 permeate had accumulated. Continue concentrating to minimum vol. Stop at 17:37. Remove 3.266g (~3.15 mL) retentate. Rinse apparatus with 6 mL water, save separately.

Begin stirring retentate in ice bath. Remove 110  $\mu$ l sample + dilute with 140  $\mu$ l (1/20) carbonate HPLC buffer at ~10:15. NPC-10K-PEG weighed yesterday (344 mg initially, dropping to 337 mg over 1-2 min static electricity?) + stored covered in freezer. Added 337 mg x 4  $\mu$ l x 90% = 1.22 mL 1 mM HCE at 10:15. Again, yellow color on solid until well suspended. Refrigerate until add to rxn dropwise (~1 drop/sec) 10:59:23 - 11:04:10. Rinse tube with 0.12 mL mM HCE after reducing foam under partial vacuum. Add rxn to rxn at 11:07. Add 2 mL carbonate HPLC buffer to PEG-tube.  $UV_{400 \text{ nm}} @ 13:38 = \frac{1.01}{18.4 \text{ g/sec/min}} = 0.055 \text{ mm} \times 10 \frac{\text{mg}}{\text{mL}}$   
 At 11:18:00 remove 20  $\mu$ l rxn sample for W monitoring at room temp.  $t_{1/2} = 11.5 \text{ min}$  (see LDW 01214.KD). Drop 0.1 mm W cell, breaking it. Dilute  $\mu$ l V soln all recovered (<20  $\mu$ l) V soln with 380  $\mu$ l water (~1/20) for HPLC 9901214.01. PEG submit calcd as 7.5, so

Continued on Page 5

Read and Understood By

L D Williams

Signed

Date

Mark Sifer

Signed

Date

PROJECT \_\_\_\_\_

quench 2 mL in 22 mL 260 mM Gly Gly ( )  
Filter through 0.45  $\mu$  filter. See [redacted] & OI for HPLC.

~~Start ultrafiltration on 30K Membrane~~

Took wrong tube from refrigerator: ultrafiltering Round 1 permeate  
No activity in retentate or permeate. Ultrafiltered 0.1 mg/mL BSA/50 mM  
borax to try to chase off bound activity before discovering there was no  
activity to chase off.

[redacted] Round 3: Removed 10  $\mu$ L sample of unquenched round 2 rxn  
for diln with 190  $\mu$ L water + HPLC.

Dissolved 102.4 mg NPC-10K-PEG in 0.36 mL 1 mM HCl and added  
to the round 2 rxn at 17:23 with stirring in an ice bath. Rinsed in  
PEG with 40  $\mu$ L 1 mM HCl.

EXHIBIT A PAGE 4 OF 4

Continued on Page

Read and Understood By

*Mark Sawyer*

Signed

Date

Signed

Date